

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Mouse Alzheimer's Disease

Cat. no. 330231 PAMM-057ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® Mastercycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000®
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche® LightCycler® 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm® BioMark™



Sample & Assay Technologies

Description

The Mouse Alzheimer's Disease RT² Profiler™ PCR Array profiles the expression of 84 genes important in the onset, development, and progression of Alzheimer's disease. The array includes genes that contribute to amyloid beta-peptide (Aβ) generation, clearance, and degradation, as well as genes involved in amyloid beta-peptide (Aβ) signal transduction leading to neuronal toxicity and inflammation. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to Alzheimer's Disease with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C.

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	A2m	Abca1	Ache	Adam9	Apba1	Apba3	Apbb1	Apbb2	Aph1a	Aplp1	Aplp2	Apoa1
B	Apoe	App	Bace1	Bace2	Bche	Bdnf	Casp3	Casp4	Cdk1	Cdk5	Cdkl1	Chat
C	Clu	Ctsb	Ctsc	Ctsd	Ctsg	Ctsl	Ep300	Ern1	Gap43	Gnao1	Gnaz	Gnb1
D	Gnb2	Gnb4	Gnb5	Gng10	Gng11	Gng3	Gng4	Gng5	Gng7	Gng8	Gngt1	Gngt2
E	Gsk3a	Gsk3b	Hsd17b10	Ide	Igf2	Ii1a	Insr	Lpl	Lrp1	Lrp6	Lrp8	Mapt
F	Mpo	Mtap2	Nae1	Ncstn	Pkp4	Plat	Plau	Plg	Prkca	Prkcb	Prkcc	Prkcd
G	Prkce	Prkci	Prkq	Prkcz	Psen1	Psen2	Serpina3c	Snca	Sncb	Ubqln1	Uqcrc1	Uqcrc2
H	Actb	B2m	Gapdh	Gusb	Hsp90ab1	MGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mm.30151	NM_175628	A2m	Alpha-2-macroglobulin
A02	Mm.277376	NM_013454	Abca1	ATP-binding cassette, sub-family A (ABC1), member 1
A03	Mm.255464	NM_009599	Ache	Acetylcholinesterase
A04	Mm.28908	NM_007404	Adam9	A disintegrin and metallopeptidase domain 9 (meltrin gamma)
A05	Mm.22879	NM_177034	Apba1	Amyloid beta (A4) precursor protein binding, family A, member 1
A06	Mm.293931	NM_018758	Apba3	Amyloid beta (A4) precursor protein-binding, family A, member 3
A07	Mm.38469	NM_009685	Apbb1	Amyloid beta (A4) precursor protein-binding, family B, member 1
A08	Mm.5159	NM_009686	Apbb2	Amyloid beta (A4) precursor protein-binding, family B, member 2
A09	Mm.359706	NM_146104	Aph1a	Anterior pharynx defective 1a homolog (C. elegans)
A10	Mm.2381	NM_007467	Aplp1	Amyloid beta (A4) precursor-like protein 1
A11	Mm.19133	NM_009691	Aplp2	Amyloid beta (A4) precursor-like protein 2
A12	Mm.26743	NM_009692	Apoa1	Apolipoprotein A-I
B01	Mm.305152	NM_009696	Apoe	Apolipoprotein E
B02	Mm.277585	NM_007471	App	Amyloid beta (A4) precursor protein
B03	Mm.24044	NM_011792	Bace1	Beta-site APP cleaving enzyme 1
B04	Mm.97885	NM_019517	Bace2	Beta-site APP-cleaving enzyme 2
B05	Mm.250719	NM_009738	Bche	Butyrylcholinesterase
B06	Mm.1442	NM_007540	Bdnf	Brain derived neurotrophic factor
B07	Mm.34405	NM_009810	Casp3	Caspase 3
B08	Mm.1569	NM_007609	Casp4	Caspase 4, apoptosis-related cysteine peptidase
B09	Mm.281367	NM_007659	Cdk1	Cyclin-dependent kinase 1
B10	Mm.298798	NM_007668	Cdk5	Cyclin-dependent kinase 5
B11	Mm.132325	NM_183294	Cdkl1	Cyclin-dependent kinase-like 1 (CDC2-related kinase)
B12	Mm.442817	NM_009891	Chat	Choline acetyltransferase
C01	Mm.200608	NM_013492	Clu	Clusterin
C02	Mm.236553	NM_007798	Ctsb	Cathepsin B
C03	Mm.322945	NM_009982	Ctsc	Cathepsin C
C04	Mm.231395	NM_009983	Ctsd	Cathepsin D
C05	Mm.4858	NM_007800	Ctsg	Cathepsin G
C06	Mm.930	NM_009984	Ctsl	Cathepsin L
C07	Mm.258397	NM_177821	Ep300	E1A binding protein p300
C08	Mm.340943	NM_023913	Ern1	Endoplasmic reticulum (ER) to nucleus signalling 1
C09	Mm.1222	NM_008083	Gap43	Growth associated protein 43
C10	Mm.251445	NM_010308	Gnao1	Guanine nucleotide binding protein, alpha o
C11	Mm.393508	NM_010311	Gnaz	Guanine nucleotide binding protein, alpha z subunit
C12	Mm.2344	NM_008142	Gnb1	Guanine nucleotide binding protein (G protein), beta 1
D01	Mm.30141	NM_010312	Gnb2	Guanine nucleotide binding protein (G protein), beta 2
D02	Mm.139192	NM_013531	Gnb4	Guanine nucleotide binding protein (G protein), beta 4
D03	Mm.17604	NM_010313	Gnb5	Guanine nucleotide binding protein (G protein), beta 5
D04	Mm.41780	NM_025277	Gng10	Guanine nucleotide binding protein (G protein), gamma 10
D05	Mm.25547	NM_025331	Gng11	Guanine nucleotide binding protein (G protein), gamma 11
D06	Mm.482164	NM_010316	Gng3	Guanine nucleotide binding protein (G protein), gamma 3
D07	Mm.215394	NM_010317	Gng4	Guanine nucleotide binding protein (G protein), gamma 4
D08	Mm.140804	NM_010318	Gng5	Guanine nucleotide binding protein (G protein), gamma 5
D09	Mm.222496	NM_010319	Gng7	Guanine nucleotide binding protein (G protein), gamma 7

Position	UniGene	GenBank	Symbol	Description
D10	Mm.103770	NM_010320	Gng8	Guanine nucleotide binding protein (G protein), gamma 8
D11	Mm.95398	NM_010314	Gngt1	Guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 1
D12	Mm.46299	NM_023121	Gngt2	Guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 2
E01	Mm.476745	NM_001031667	Gsk3a	Glycogen synthase kinase 3 alpha
E02	Mm.394930	NM_019827	Gsk3b	Glycogen synthase kinase 3 beta
E03	Mm.6994	NM_016763	Hsd17b10	Hydroxysteroid (17-beta) dehydrogenase 10
E04	Mm.28366	NM_031156	Ide	Insulin degrading enzyme
E05	Mm.3862	NM_010514	Igf2	Insulin-like growth factor 2
E06	Mm.15534	NM_010554	Il1a	Interleukin 1 alpha
E07	Mm.268003	NM_010568	Insr	Insulin receptor
E08	Mm.1514	NM_008509	Lpl	Lipoprotein lipase
E09	Mm.271854	NM_008512	Lrp1	Low density lipoprotein receptor-related protein 1
E10	Mm.321990	NM_008514	Lrp6	Low density lipoprotein receptor-related protein 6
E11	Mm.442134	NM_001080926	Lrp8	Low density lipoprotein receptor-related protein 8, apolipoprotein e receptor
E12	Mm.1287	NM_010838	Mapt	Microtubule-associated protein tau
F01	Mm.4668	NM_010824	Mpo	Myeloperoxidase
F02	Mm.256966	NM_001039934	Mtap2	Microtubule-associated protein 2
F03	Mm.237670	NM_144931	Nae1	NEDD8 activating enzyme E1 subunit 1
F04	Mm.218203	NM_021607	Ncstn	Nicastrin
F05	Mm.260938	NM_026361	Pkp4	Plakophilin 4
F06	Mm.154660	NM_008872	Plat	Plasminogen activator, tissue
F07	Mm.4183	NM_008873	Plau	Plasminogen activator, urokinase
F08	Mm.971	NM_008877	Plg	Plasminogen
F09	Mm.222178	NM_011101	Prkca	Protein kinase C, alpha
F10	Mm.207496	NM_008855	Prkcb	Protein kinase C, beta
F11	Mm.7980	NM_011102	Prkcc	Protein kinase C, gamma
F12	Mm.2314	NM_011103	Prkcd	Protein kinase C, delta
G01	Mm.24614	NM_011104	Prkce	Protein kinase C, epsilon
G02	Mm.291554	NM_008857	Prkci	Protein kinase C, iota
G03	Mm.329993	NM_008859	Prkck	Protein kinase C, theta
G04	Mm.28561	NM_008860	Prkcz	Protein kinase C, zeta
G05	Mm.998	NM_008943	Psen1	Presenilin 1
G06	Mm.330850	NM_011183	Psen2	Presenilin 2
G07	Mm.422814	NM_008458	Serpina3c	Serine (or cysteine) peptidase inhibitor, clade A, member 3C
G08	Mm.17484	NM_009221	Snca	Synuclein, alpha
G09	Mm.200843	NM_033610	Sncb	Synuclein, beta
G10	Mm.182053	NM_026842	Ubqln1	Ubiquilin 1
G11	Mm.335460	NM_025407	Uqcrc1	Ubiquinol-cytochrome c reductase core protein 1
G12	Mm.334206	NM_025899	Uqcrc2	Ubiquinol cytochrome c reductase core protein 2
H01	Mm.328431	NM_007393	Actb	Actin, beta
H02	Mm.163	NM_009735	B2m	Beta-2 microglobulin
H03	Mm.343110	NM_008084	Gapdh	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mm.3317	NM_010368	Gusb	Glucuronidase, beta
H05	Mm.2180	NM_008302	Hsp90ab1	Heat shock protein 90 alpha (cytosolic), class B member 1
H06	N/A	SA_00106	MGDC	Mouse Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies